




INSTRUCTIONS FOR USE FOR THE SOFT TISSUE

CE 0426

ASPIRATION CUTTING BIOPSY INSTRUMENT

ATA-CUT[®]

(GELMAN NEEDLE OR MODIFIED MENGHINI)



- WARNING:**
- Store in a dry clean place away from heat sources
 - Before use ensure that the packaging is intact: if perforated, damaged or open discard the instrument
 - Check the expiry date and that the model corresponds
 - Peel open the package
 - Pay attention to not compromise the sterility of the instrument before use
 - In the event of a product defect or malfunction complete the yellow claim form to be found in the product carton and conserve the defective instrument

N.B.:

The use of this instrument is reserved to qualified medical personnel who are aware of its application.

The instrument is sterile (sterilized by ETO GAS) single use and should be discarded after use (Circ.Min:Sanità 403/9.5/1146-L. 10.81)

Sterylab declines any and all responsibility in the event of improper use or reutilisation.

BIOPTIC PROCEDURE

- 1) Push the piston to the bottom of the syringe depressing the tongue of the self-locking system
- 2) Pencil-grip the instrument and reach the lesion (organ) to biopsy.

In the event of an oblique penetration turn so that the self-locking system is pointing upwards (the Internal Echogenic Marker will facilitate the visualization of the needle on the monitor).
- 3) Having reached the lesion create a depression in the syringe by pulling the piston back until the self-locking system engages
- 4) Penetrate into the lesion and section the tissue
- 5) Remove the needle from the patient
- 6) Expel the histological cylinder by pushing the piston to the bottom of the syringe after depressing the tongue of the self-locking system