



## Seed Alignment Tray





## **Table of Contents**


General Precautions .....	2
Customer Responsibility .....	3
Operation .....	4
Maintenance .....	5
Service .....	5
Parts and Accessories List .....	5
Features and Specifications .....	5
Warranty .....	6

## **General Precautions**

 **CAUTION:** Proper use of this device depends on careful reading of all instructions and labels.

 **CAUTION:** Steam Sterilize (autoclave) or ETO sterilize only. Do not gamma ray sterilize.

 **WARNING:** After steam sterilization, unit will be hot. Allow adequate cooling time before handling. Do not drop or mishandle unit.

 **WARNING:** For use with gamma and x-ray emitting LDR radiation sources only.

## **Customer Responsibility**

⚠ Federal law in the U.S.A. and Canada restricts the sale, distribution or use of this device to, by or on the order of a licensed medical practitioner. The use of this device should be restricted to the supervision of a qualified medical physicist. Handling of radioactive sources is potentially hazardous and should be performed by qualified personnel. Should repair or replacement of this device become necessary after the warranty period, the customer should seek advice from Standard Imaging Inc. prior to such repair or replacement. If this device is in need of repair, it should not be used until all repairs have been made and the product is functioning properly and ready for use. The owner of this device has sole responsibility for any malfunction resulting from abuse, improper use or maintenance, or repair by anyone other than Standard Imaging Inc. The information in this manual is subject to change without notice. No part of this manual may be copied or reproduced in any form, or by any means, without prior written consent of Standard Imaging Inc.

## Operation

The Seed Alignment Tray is used with a SeedVac™ to conveniently load needles and reduce radiation exposure. The Seed Alignment Tray can be used alone or connected to the Seed Sterilization and Sorting Tray.

1. Attach the Seed Alignment Tray to the Loading Cover of the Seed Sterilization and Sorting Tray by removing the two screws on top of the Loading Cover. Align the captive screws on the Alignment Tray to the nuts revealed in the Loading Cover. Secure finger-tight.



*Fig. 1: The Seed Alignment Tray attached to the Loading Cover.*

2. Sterilize the assembly with the loaded Sterilization Tray and spatula.
3. Remove the Sterilization Cover and attach the Loading Cover.
4. Use the sterilized spatula to scoop out the seeds from the wells.
5. Pour spacers onto the Seed Alignment Tray.

6. Align the seeds and spacers in the slots of the tray.
7. Use a SeedVac to pick up the seeds and spacers.
8. Proceed with loading the needles.

## **Maintenance**

As is standard practice, it is recommended that the Seed Alignment Tray be examined after each use for deterioration. To remove any dirt, use a cloth slightly dampened with a solution of mild detergent and water or a liquid disinfecting agent.

## **Service**

There are no serviceable parts on the tray.

## **Parts and Accessories**

<b>REF</b>	<b>Description</b>
90089	Seed Alignment Tray
80201	Instruction Manual

## **Features and Specifications**

Material .....	Aluminum
Width .....	4.5 in. (11.4 cm)
Depth .....	4.5 in. (11.4 cm)

*Specifications subject to change without notice.*

## Warranty

This product is sold by Standard Imaging Inc. under the warranty herein set forth. The warranty is extended only to the buyer purchasing the product directly from Standard Imaging Inc. or as a new product from an authorized dealer or distributor of Standard Imaging Inc.

For a period of twelve (12) months from the date of original delivery to the purchaser or a distributor, this product is warranted against functional defects in materials and workmanship, provided it is properly operated under conditions of normal use, and that repairs and replacements are made in accordance herewith. The foregoing warranty shall not apply if the product has been disassembled, altered or repaired other than by Standard Imaging Inc. or if the product has been subject to abuse, misuse, negligence or accident.

Standard Imaging's sole and exclusive obligation and the purchaser's sole and exclusive remedy under the above warranties are limited to repairing or replacing free of charge, at Standard Imaging's option, a product: (1) which contains a defect covered by the above warranties; (2) which are reported to Standard Imaging not later than seven (7) days after the expiration date of the 12-month warranty period; (3) which are returned to Standard Imaging promptly after discovery of the defect; and (4) which are found to be defective upon Standard Imaging's examination. Transportation charges are the buyer's responsibility. STANDARD IMAGING INC. SHALL NOT BE OTHERWISE LIABLE FOR ANY DAMAGES, INCLUDING BUT NOT LIMITED TO, INCIDENTAL DAMAGES, CONSEQUENTIAL DAMAGES, OR SPECIAL DAMAGES.

THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL STANDARD IMAGING INC. BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM THE USE, MISUSE OR ABUSE OF THE PRODUCT OR CAUSED BY ANY DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM OF SUCH DAMAGE IS BASED UPON THE WARRANTY, CONTRACT, NEGLIGENCE, OR OTHERWISE.

Distributed by SeeDOS Ltd  
For further information or a quotation please  
contact Colin Walters at [cwalters@seedos.com](mailto:cwalters@seedos.com)