

Hawkins™ Blunt Needle Access System

Hawkins™ Blunt Needle Access System

- Provides a lightweight method of targeting, remaining in place while position is confirmed. This is especially valuable in CT guided procedures.
- The adjustable needle stop localises the outer cannula to the target, preventing forward migration.
- The outer cannula has an echogenic tip, promoting accurate placement under ultrasound guidance.
- Numerically ordered centimetre markings on the outer cannula facilitate precise depth placement.
- Accommodates most guidewires, biopsy needles, and biopsy instruments.

Hawkins™ Blunt Needle Access System

		Box of 5
218057	18G (1.2 mm) x 5 cm	ID 1.0 mm
218107	18G (1.2 mm) x 10 cm	ID 1.0 mm
218157	18G (1.2 mm) x 15 cm	ID 1.0 mm
218207	18G (1.2 mm) x 20 cm	ID 1.0 mm
216057	16G (1.6 mm) x 5 cm	ID 1.3 mm
216107	16G (1.6 mm) x 10 cm	ID 1.3 mm
216157	16G (1.6 mm) x 15 cm	ID 1.3 mm

218157

