

## PROSTATE LOADING SYSTEMS

### Seed Vac Handling and Loading Starter Kit- Model SHK.4

90091	Seed Vac
90085	Seed Sterilization and Sorting Tray
72012	Shielding Cover for Remote Sterilization with 90085
90072	Needle Loading Platform with Needle Loading Shield
90080	Case of Ten Protective Lead Seed Pouch Kits
90089	Seed Alignment Tray
96823A1	Ruler, Chrome Machinists' 15cm; 1/2mm reading front and back
90040-G	Reverse Action 6.5 inch Tweezers with Groove, Single Tweezers
99503	Poly-Vac Sterilization Tray
99017	Portable Survey Meter Model 3, suitable for LDR brachytherapy seeds includes 44-3 Gamma Scintillation Detector Meter read out is in counts per minute and microR/hr
99816	Diddler (Needle Adjustment Tool)
99053I	Tissue Equivalent Ultrasound Practice Phantom
98212	Dummy Seeds, Pack of 50
91514B LOOK	Pre-cut Spacers, 900 pieces
90065	Seed Sensor

### Seed Slider Handling and Loading Starter Kit- Model SHK.3

90090	Seed Slider
90085	Seed Sterilization and Sorting Tray
72012	Shielding Cover for Remote Sterilization with 90085
90070	Needle Loading Shield
90080	Case of Ten Protective Lead Seed Pouch Kits
90095	Needle Cradle
96823A1	Ruler, Chrome Machinists' 15cm; 1/2mm reading front and back
90040-G	Reverse Action 6.5 inch Tweezers with Groove Single Tweezers
99503	Poly-Vac Sterilization Tray
99017	Portable Survey Meter Model 3, suitable for LDR brachytherapy seeds includes 44-3 Gamma Scintillation Detector Meter read out is in counts per minute and microR/hr
98816	Diddler (Needle Adjustment Tool)
99053I	Tissue Equivalent Ultrasound Practice Phantom
98212	Dummy Seeds, Pack of 50
91514B	LOOK Pre-cut Spacers, 900 pieces
90065	Seed Sensor

### MICK Seed Handling and Loading Starter Kit- Model SHK.M

90087	MICK Magazine Sterilizer
90070	Needle Loading Shield
90091	SeedVac
90080	Case of Ten Protective Lead Seed Pouch Kits

# Brachytherapy LDR Calibration Systems

## **BRS-4A For: Verification of $^{125}\text{I}$ and $^{103}\text{Pd}$ Seed Strengths Individually and in Batches of up to 500 Seeds**

Consisting of:

- Premier 3000 Electrometer
- ADCL Electrometer Calibration in Amps and Coulombs
- HDR 1000 Plus Ion Chamber
- Source Holder for Single LDR Seed Assay, P/N 70016
- ADCL Calibration for One Iodine or Palladium Isotope
- Source Holder for Batch LDR Seed Assay, P/N 70022
- Carrying Case

## **BRS-5A For: Verification of I-125 RAPID Strand™**

Consisting of:

- MAX 4000 Electrometer
- ADCL Calibration for Electrometer
- HDR 1000 Plus Ion Chamber
- Source Holder for Single LDR Seed Assay, P/N 70016
- ADCL Calibration I-125 RAPID Strand™ code IMC6711, P/N 80040-B
- Source Holder for I-125 RAPID Strand™ Iodine Seeds, P/N 70023
- Carrying Case

## **BRS-6A For: Verification of Indigo™ EXPRESS Seeding Cartridge**

Consisting of:

- MAX 4000 Electrometer
- ADCL Calibration for Electrometer
- HDR 1000 Plus Ion Chamber
- Source Holder for Single LDR Seed Assay, P/N 70016
- ADCL Calibration for Theragenics Theraseed 200, P/N 80035-B
- Source Holder for Indigo™ EXPRESS Seeding Cartridge, P/N 70032
- Carrying Case

Distributed by SeedOS Ltd

For further information and a quotation please contact  
Colin Walters at [cwalters@seedos.com](mailto:cwalters@seedos.com)  
[www.seedos.com](http://www.seedos.com)

Doc No. 1129-03  
10-31-01, page 2 of 2